

Serum Iron Studies

CPT Code 82728, 83540, 83550, 84466



CMS Policy

Local policies are determined by the performing test location. This is determined by the state in which your performing laboratory resides and where your testing is commonly performed.

Medically Supportive ICD Codes are listed on subsequent page(s) of this document.

Coverage Indications, Limitations, and/or Medical Necessity

Serum iron studies are effective in the evaluation of disorders of iron metabolism, iron deficiency and iron excess. For best results, iron studies should be performed when the patient is fasting in the morning and has abstained from medications that may influence iron balance. The most common cause of anemia is iron deficiency. In children on a milk diet, iron deficiency is typically secondary to dietary deficiency. In adults, iron deficiency is usually the result of blood loss and is only occasionally secondary to dietary deficiency or malabsorption. Patient's may have iron deficient erythropoietin for months or years if adequate iron replacement has not been given following surgery. Serum iron may be elevated with high doses of supplemental iron. Serum iron can fluctuate in acute and chronic inflammatory and neoplastic conditions.

Total Iron Binding Capacity (TIBC) is an indirect measure of transferrin, a protein that binds and transports iron. TIBC quantifies transferrin by the amount of iron that it can bind. TIBC and transferrin are elevated in iron deficiency, and with oral contraceptive use, and during pregnancy. TIBC and transferrin may be decreased in malabsorption syndromes or in those affected with chronic diseases. The percent saturation represents the ratio of iron to the TIBC.

Assays for ferreting are also useful in assessing iron balance. Low concentrations are associated with iron deficiency and are highly specific. High concentrations are found in hemosiderosis (iron overload without associated tissue injury) and hemochromatosis (iron overload with associated tissue injury). In these conditions the iron is elevated, the TIBC and transferrin are within the reference range or low, and the percent saturation is elevated. Serum ferritin can be useful for both initiating and monitoring treatment for iron overload.

Transferrin and ferritin belong to a group of serum proteins known as acute phase reactants, and are increased in response to stressful or inflammatory conditions and also can occur with infection and tissue injury due to surgery, trauma or necrosis. Ferritin and iron/TIBC (or transferrin) are affected by acute and chronic inflammatory conditions, and in patients with these disorders, tests of iron status may be difficult to interpret.

Indications

1. Ferritin, iron and either iron binding capacity or transferrin are useful in the differential diagnosis of iron deficiency, anemia, and for iron overload conditions.

- a. The following presentations are examples that may support the use of these studies for evaluating iron deficiency:
- ◆ Certain abnormal blood count values (i.e., decreased Mean Corpuscular Volume (MCV), decreased hemoglobin/hematocrit when the MCV is low or normal, or increased red cell Distribution Width (RDW) and low or normal MCV) Abnormal appetite (pica)
 - ◆ Acute or chronic gastrointestinal blood loss
 - ◆ Hematuria
 - ◆ Menorrhagia
 - ◆ Malabsorption
 - ◆ Status post-gastrectomy
 - ◆ Status post-gastrojejunostomy
 - ◆ Malnutrition
 - ◆ Preoperative autologous blood collection(s)
 - ◆ Malignant, chronic inflammatory and infectious conditions associated with anemia which may present in a similar manner to iron deficiency anemia
 - ◆ Following a significant surgical procedure where blood loss had occurred and had not been repaired with adequate iron replacement

Medicare Local Coverage Determination Policy

Serum Iron Studies

CPT Code 82728, 83540, 83550, 84466



- b. The following presentations are examples that may support the use of these studies for evaluating iron overload:
- ◆ Chronic Hepatitis
 - ◆ Diabetes
 - ◆ Hyperpigmentation of skin
 - ◆ Arthropathy
 - ◆ Cirrhosis
 - ◆ Hypogonadism
 - ◆ Hypopituitarism
 - ◆ Impaired porphyrin metabolism
 - ◆ Heart failure
 - ◆ Multiple transfusions
 - ◆ Sideroblastic anemia
 - ◆ Thalassemia major
 - ◆ Cardiomyopathy, cardiac dysrhythmias and conduction disturbances
2. Follow-up testing may be appropriate to monitor response to therapy, e.g., oral or parenteral iron, ascorbic acid, and erythropoietin.
3. Iron studies may be appropriate in patients after treatment for other nutritional deficiency anemias, such as folate and vitamin B12, because iron deficiency may not be revealed until such a nutritional deficiency is treated.
4. Serum ferritin may be appropriate for monitoring iron status in patients with chronic renal disease with or without dialysis.
5. Serum iron may also be indicated for evaluation of toxic effects of iron and other metals (e.g., nickel, cadmium, aluminum, and lead) whether due to accidental, intentional exposure or metabolic cause.

Limitations

1. Iron studies should be used to diagnose and manage iron deficiency or iron overload states. These tests are not to be used solely to assess acute phase reactants where disease management will be unchanged. For example, infections and malignancies are associated with elevations in acute phase reactants such as ferritin, and decreases in serum iron concentration, but iron studies would only be medically necessary if results of iron studies might alter the management of the primary diagnosis or might warrant direct treatment of an iron disorder or condition.
2. If a normal serum ferritin level is documented, repeat testing would not ordinarily be medically necessary unless there is a change in the patient's condition, and ferritin assessment is needed for the ongoing management of the patient. For example, a patient presents with new onset insulin-dependent diabetes mellitus and has a serum ferritin level performed for the suspicion of hemochromatosis. If the ferritin level is normal, the repeat ferritin for diabetes mellitus would not be medically necessary.
3. When an End Stage Renal Disease (ESRD) patient is tested for ferritin, testing more frequently than every three months requires documentation of medical necessity (e.g., other than chronic renal failure or renal failure, unspecified).
4. It is ordinarily not necessary to measure both transferrin and TIBC at the same time because TIBC is an indirect measure of transferrin. When transferrin is ordered as part of the nutritional assessment for evaluating malnutrition, it is not necessary to order other iron studies unless iron deficiency or iron overload is suspected as well.
5. It is not necessary to measure either iron/TIBC (or transferrin) and ferritin in initial patient testing. If clinically indicated after evaluation of the initial iron studies, it may be appropriate to perform additional iron studies either on the initial specimen or on a subsequently obtained specimen. After a diagnosis of iron deficiency or iron overload is established, either iron/TIBC (or transferrin) or ferritin may be medically necessary for monitoring, but not both.
6. It would not ordinarily be considered medically necessary to do a ferritin as a preoperative test except in the presence of anemia or recent autologous blood collections prior to the surgery
-

Medicare Local Coverage Determination Policy

Serum Iron Studies

CPT Code 82728, 83540, 83550, 84466



The ICD10 codes listed below are the top diagnosis codes currently utilized by ordering physicians for the limited coverage test highlighted above that are also listed as medically supportive under Medicare's limited coverage policy.

***Note—Bolded diagnoses below have the highest utilization**

Code	Description
D50.0	Iron deficiency anemia secondary to blood loss (chronic)
D50.8	Other iron deficiency anemias
D50.9	Iron deficiency anemia, unspecified
D51.0	Vitamin B12 deficiency anemia due to intrinsic factor deficiency
D51.8	Other vitamin B12 deficiency anemias
D51.9	Vitamin B12 deficiency anemia, unspecified
D52.9	Folate deficiency anemia, unspecified
D53.9	Nutritional anemia, unspecified
D63.1	Anemia in chronic kidney disease
D63.8	Anemia in other chronic diseases classified elsewhere
D64.9	Anemia, unspecified
E11.22	Type 2 diabetes mellitus with diabetic chronic kidney disease
E11.65	Type 2 diabetes mellitus with hyperglycemia
E11.9	Type 2 diabetes mellitus without complications
E61.1	Iron deficiency
M25.50	Pain in unspecified joint
N18.4	Chronic kidney disease, stage 4 (severe)
N18.9	Chronic kidney disease, unspecified
R79.89	Other specified abnormal findings of blood chemistry
R79.9	Abnormal finding of blood chemistry, unspecified

Disclaimer:

This diagnosis code reference guide is provided as an aid to physicians and office staff. Diagnosis codes must be applicable to the patient's symptoms or conditions and must be consistent with documentation in the patient's medical record. Kashi Clinical Laboratories does not recommend any diagnosis codes and will only submit diagnosis information provided to us by the ordering physician or his/her designated staff. The CPT codes provided are based on AMA guidelines and are for informational purposes only. CPT coding is the sole responsibility of the billing party. Please direct any questions regarding coding to the payer being billed.